

WS-YDB Horizontal Cylindrical Pressure Steam Sterilizer



Characteristics

1. Sterilizing course automatic sterilization controlled by computer, easy to operate.
2. With drying function, suitable dressing drying.
3. With over-temperature, over-pressure auto-protect device.
4. The door opening mechanism cannot be operated until the pressure in the chamber is reduced to 0.027MPa. It can't be start on if the doesn't closed well.
5. The safety value will be open automatically when the inner pressure over 0.24MPa, and the steam be exhausting to the water tank.
6. The power will be cut-off automatically, cut-off the water and alarm when the machine lacking water.
7. The surface of the sterilizer is steel coated, the chamber of the sterilizer is made of stainless steel.
8. With printer.

Technical Data

Model technical data		WS-280YDB	WS-200YDB	WS-150YDB
Sterilizing	chamber	280L($\phi 600 \times 1000$ mm)	200L($\phi 515 \times 1000$ mm)	150L($\phi 440 \times 1000$ mm)

volume			
working pressure	0.22 MPa		
working temperature	134℃		
Adjustment of temperature	40℃ - 134℃		
Time of sterilization	0-99 min		
Time of drying	0-99 min		
Heat average	≤± 2℃		
Power	12KW/380V 50Hz	9KW/380V 50Hz	9KW/380V 50Hz
Dimension	1520×910×1900	1400×750×1700	1400×650×1600
Package	1680×1080×2100	1560×910×1880	1560×820×1820
GW/NW	570/462 KG	436/350 KG	430/340 KG