

DBP-3000C VCM Doppler Blood Pressure Monitor



Introduction

Doppler Blood Pressure is used to test the animal blood pressure in hospital.

Features

- Aneroid Sphygmomanometer, no pollution
- Integrated sphygmomanometer, easy to operate
- Specialized probe with high sensitivity
- Abundant vascular doppler menu
- Multiple size cuff, feasible to different animals
- Display vascular sound volume and parameter
- Turn off when no using for 3 minutes
- Both of battery and AC power supply

Technical Specifications

PARAMETER	SPECIFICATION
Display	Color display,support touch screen
Blood pressure range	0-300mmHg
Pressure precision	±3.75mmHg
Vascular probe frequency	8Mhz
Powe supply	Bulit-in rechargeable battery AC220/110V 50/60HZ

Configuration

Website: www.hqs-vetcare.com

Add: Room 402, Unit 6, Building 27, Pingchou District 2, Jiangdong Street,
Yiwu, Jinhua, Zhejiang Province

1pc Vascular doppler
1pc 5/8Mhz vascular probe
1pc Popular Sphygmomanometer
1set Cuff
1pc Adaptor
1pc Carry case
Optional Welchallyn sphygmomanometer